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ACERO MONOFILAMIENTO SIN AGUJA: 10372360019

ACERO MONOFILAMIENTO CON AGUJA: 10372360007

INSTRUCTIONS FOR USE OF MONOFILAMENT STEEL

NON ABSORBABLE SURGICAL SUTURE

DESCRIPTION

Monofilament Steel Suture is a monofilament, synthetic, non absorbable surgical suture, composed of 316L stainless steel. Monofilament Steel Suture meets such specifications as set forth for Non Absorbable Surgical Sutures by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia and for the ABNT norm NBR 13904/2003 - Filaments for Surgical Sutures. It is sterilized by Gamma Radiations (Cobalt 60) or Ethylene Oxide. Monofilament Steel Suture meet all requirements established by the norms for non absorbable surgical sutures, except the following diameter:

MÁXIMUM SUTURE OVERSIZE IN DIAMETER (mm) FROM U.S.P	
U.S.P - SUTURE SIZE DESIGNATION	MAXIMUM OVERSIZE (mm)
5	0.001

INDICATION

Monofilament Steel Suture is indicated for use in abdominal wound closure, hernia repair, sternal closure and orthopedic procedures, including cases of bone and restoration of tendons. Monofilament Steel Suture Without Needle is indicated for procedures in maxillo-facial and can be used in the procedure of neuro surgery and orthopedics.

ACTION

Monofilament Steel Suture elicits a minimal acute inflammatory reaction in tissue and is not absorbed.

CONTRAINDICATION

Monofilament Steel Suture is contraindicated in patients with known sensitivities or allergies to 316L stainless steel, or constituent metals such as chromium and nickel.

WARNINGS

Users should be familiar with surgical procedures and techniques involving non absorbable sutures before employing Monofilament Steel Suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Acceptable surgical practice should be followed for the management under surgeon's orientation according to applicable surgical techniques and it should only be manipulated by qualified people. Do not use the suture for other purposes, just the one indicated by the manufacturer. When handling this suture be sure of doing it aseptically. Acceptable surgical practice should be followed for the management of infected or contaminated wounds. Do not re-sterilize. Discard the open and unused packages. Discard left over. Do not use if packing is violated or damp. Do not use the product beyond the validity period. Disposal of used sutures: Risk of biological contamination - adopt internal disposal procedures, taking into account the risk of biological contamination. Disposal of used and non-used sutures: Perforating and cutting materials - adopt internal disposal procedures for perforating and cutting materials.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Braking of needles may cause surgeries to take a longer, additional surgeries and foreign matters left inside the body. Inadvertent piercing of users with contaminated needles may cause transmission of pathogens through blood. Users must be careful when handling surgical needles to avoid undesirable perforations. Dispose used needles in recipients for potentially contaminated material.

USE

Use the necessary amount of suture in agreement with the extension of the surgical wound and the technique adopted by the surgeon. The size and diameter of suture and the needle used for each type of surgical intervention are defined as per medical criteria.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, gradual loss of tensile strength over time, allergic response in patients with know sensitivities to 316L stainless steel, or constituent metals such as chromium and nickel, infection, minimal acute inflammatory tissue reaction, pain, edema and local irritation at the wound site.

STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30°C (59°F to 86°F), protected from direct sunlight, heat sources, humidity and dust.

HOW SUPPLIED

Monofilament Steel Suture is supplied in natural color, in sterile envelopes in U.S.P sizes 7-0 through 5 (metric 0.7 through 9). Suture is supplied in pre-established lengths, without needle or attached with 01 (one) or 02 (two) non traumatic needle, according to each type of surgical intervention; in boxes with 12, 24 or 36 units of units.

EXPIRATION

Five years from manufacturing date. See box and envelopes.

PACKAGES WITH 12 (TWELVE), 24 (TWENTY FOUR) OR

36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.

SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES

DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN - DESTROY AFTER USING

STERILIZATION METHOD:

GAMMA RADIATION (COBALT 60) OR ETHYLENE OXIDE

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CNPJ 02.370.649/0001-20 - Fábrica: Tel.: 55 19 3646-7300

Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

REGISTRO ANVISA N°:

MONOFILAMENT STEEL NO NEEDLE: 10372360019

MONOFILAMENT STEEL WITH NEEDLE: 10372360007