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INSTRUCTIONS FOR USE OF CHROMIC CATGUT

ABSORBABLE SURGICAL SUTURE

DESCRIPTION

Chromic Catgut Suture is an absorbable, sterile surgical suture, composed of purified connective tissue (mostly collagen) derived from the serosal layer of beef (bovine). Chromic Catgut Suture are processed to provide greater resistance to absorption and to this it is tanned by treatment with chromium salts and pyragolol, to prolong its resistance to absorption and to colour the material dark brown. Chromic Catgut Suture is packed in isopropyl alcohol solution, composed of 89.6% isopropyl alcohol p.a., 1.0% of diethanolamine, 0.25 % sodium benzoate (USP pharmaceutical grade), and sterile water injection q.s. for 100%, also called conservative solution, to maintain its pliancy. This product meets specifications set forth for Absorbable Surgical Sutures by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia (EP) and for the ABNT norm NBR13904/2003 - Filaments for Surgical Sutures. It is sterilized by Gamma Radiation (Cobalt 60).

INDICATION

Chromic Catgut suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmology, but not indicated for use in cardiovascular and neurological tissues.

ACTIONS

When Chromic Catgut suture is placed in tissue, a moderate tissue inflammation occurs which is characteristic of foreign body response to a substance. This is followed by a loss of tensile strength and a loss of suture mass, as the proteolytic enzymatic digestive process dissolves the Chromic Catgut suture. This process continues until the suture is completely absorbed. Many variable factors may affect the rate of absorption. Some of the major factors which can affect tensile strength loss and absorption rates are: type of suture (Plain Catgut suture generally absorbs more rapidly than Chromic Catgut suture); infection (Chromic Catgut suture is absorbed more rapidly infected tissue than in non-infected tissue) e tissue sites (Chromic Catgut suture will absorb more rapidly in tissue where increased levels of proteolytic enzymes are present, as in the secretions exhibited in the stomach, cervix and vagina).

CONTRAINDICATION

This suture, being absorbable, should not be used where extended approximation of tissue is required. The use of this suture is contraindicated in patients with known sensitivities or allergies to collagen or chromium, as Chromic Catgut suture is a collagen based material and it is treated with chromic salt solutions.

CARES

Users should be familiar with surgical procedures and techniques involving absorbable sutures before using Chromic Catgut suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing. As this is an absorbable material, the use of supplemental non absorbable sutures should be considered by the surgeon in the closure of sites which may undergo expansion, stretching, or distension or which may require additional support. As which any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculi formation. As an absorbable suture, Chromic Catgut suture may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of infected or contaminated wounds. Certain patients may be hypersensitive to collagen or chromium and might exhibit an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound suppuraton and bleeding, as well as sinus formation.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of surgeon. The surgeon should avoid unnecessary tension when running down knots, to reduce the occurrence of surface fraying and weakening of the strand. Avoid prolonged exposure to elevated temperatures. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

USE

Use as required for the surgical procedure.

The size and diameter of suture and needle used for each type of surgical intervention are defined as per medical criteria.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include wound dehiscence, variable rates of absorption over time (depending on such factors as the type of suture used, the presence of infection and the tissue site), failure to provide adequate wound support in closure of sites where expansion, stretching or distension occur, unless additional support is supplied through the use of non absorbable suture material, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from cancer, anemia, obesity, diabetes, infection or other conditions which may delay wound healing, allergic response in the patients with known sensitivities to collagen or chromium which may result in an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound suppuraton and bleeding, as well as sinus formation, infection, moderate tissue inflammatory response characteristic of foreign body response, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.

STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30° C (59°F to 86°F), protected from direct sunlight, heat sources, humidity and dust.

PACKING

Chromic Catgut Suture is supplied in sterile envelopes with conservative solution, in sizes 7-0 through 5 (metric sizes 0,7 through 9). Suture is supplied in pre-established lengths, without needle of with 1 or 2 non traumatic needles according to each type of surgical intervention, in boxes with 12, 24 or 36 units.

EXPIRATION

Five years from manufacturing date. See box and envelope.

PACKAGES WITH 12 (TWELVE), 24 (TWENTY FOUR) OR
36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.

SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES
DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN - DESTROY AFTER USING
STERILIZATION METHOD: GAMMA RADIATION (COBALT 60)

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CNPJ 02.370.649/0001-20 - Fábrica: Tel.: 55 19 3646-7300

Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

REGISTRO ANVISA N.º:
CHROMIC CATGUT NO NEEDLE: 10372360020
CHROMIC CATGUT WITH NEEDLE: 10372360005