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SURGICAL HEMOSTATIC

### DESCRIPTION

**Bone Wax** is a hemostatic sterile surgical of a solid mixture of whitened bee wax and Isopropyl Palmitate that acts as a softening agent. Bone Wax is slightly yellowish and has characteristic odor. This product meets specifications set forth by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia (EP). Bone Wax is sterilized by Gamma Radiations (Cobalt 60) or Ethylene Oxide.

### INDICATION

Bone wax is indicated to control bleeding on bone surfaces, such as; orthopedic surgeries, open chest surgeries where there is a cut in the sternum bone, the use is indicated for any surgeries that involve the need to cut bones.

### ACTION

Bone Wax provides local bone hemostasis, acts as mechanic barrier (plug) with no biochemical action and is slightly absorbed. It is applied in the bone bleeding region when to stop the bone bleeding is intended, till the hemostasis occurs.

### CONTRAINDICATION

Use Bone Wax is contraindicated when fast bone recovery and fusion are required.

This device is contraindicated to persons that could have any allergic reaction to the medical device material.

### WARNINGS

The users should be aware that Bone Wax may inhibit the osteogenesis and act as a physical barrier to the recovery process. Device should only be used by or under surgeon's orientation according to applicable surgical techniques and it should only be manipulated by qualified people. Do not use the Bone Wax for other purposes but those indicated by the manufacturer. Do not re-sterilize. Do not use if packing is violated or damp. Do not expose the Bone Wax to high temperatures. Do not use the product beyond the period of validity.

Disposal of used devices: Risk of biological contamination - adopt internal disposal procedures, taking into account the risk of biological contamination. Disposal of non used device can be destined to the recycling or common garbage (paper/plastic: packing) or only as common garbage (bone wax - biodegradable).

### PRECAUTIONS

Bone Wax should be used strictly and excess wax should be removed from surgical field. It should be opened immediately before the use to minimize the chances of contamination and excessive dryness.

### ADVERSE REACTIONS

Adverse effects associated with the use of Bone Wax include minimal inflammatory tissue reactions in underlying sites.

The surgeon must inform the patient to look for medical assistance in case of any suspect of adverse reactions symptoms.

### USE

Use as required for the surgical procedure. Bone Wax should be used immediately after removing it from packaging using aseptic techniques; it should be warmed by handling fingers till reaching the desired consistence.

### STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30°C (59°F to 86°F), protected from direct sunlight, heat sources, dust and humidity.

### HOW SUPPLIED

Bone Wax is available in sticks or tablets of 2.5g each, in individual sterile envelopes, in boxes with 12, 24 or 36 units.

### EXPIRATION

Five years from manufacturing date. See box and envelope.

PACKAGES WITH 12 (TWELVE), 24 (TWENTY FOUR) OR

36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.

SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES

DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN

DESTROY AFTER USING

STERILIZATION METHOD: GAMMA RADIATION OR ETHYLENE OXIDE

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CNPJ 02.370.649/0001-20 - Fábrica: Tel.: 55 19 3646-7300

Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

ANVISA REGISTRATION N.º BONE WAX: 10372360003