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NYLON MONOFILAMIENTO SIN AGUJA: 10372360014

NYLON MONOFILAMIENTO CON AGUJA: 10372360009

INSTRUCTIONS FOR USE OF NYLON MONOFILAMENT SUTURE

NON ABSORBABLE SURGICAL SUTURE

DESCRIPTION

Nylon Monofilament Suture is a synthetic, monofilament, a non absorbable, sterile, surgical suture composed of the aliphatic polymers of polyamide 6 and or polyamide 6.6. Nylon Monofilament Suture is available in black or blue to enhance its visibility in the tissue or undyed (colorless). The coloring pigments are approved by the FDA being the D&C Blue # 6 for blue nylon and Logwood Black for black nylon. The internal and external raw material suppliers may use other pigments, but they must prove that the pigment is biologically compatible and do not present toxicity. The Nylon Monofilament Suture meets such specifications as set forth for Non Absorbable Surgical Sutures by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia (EP) and for the ABNT norm NBR 13904 - Filaments for Surgical Sutures. It can be sterilized by Gamma Radiations or Ethylene Oxide.

INDICATION

The Nylon Monofilament Suture is indicated for use in the approximation and /or ligation of general soft tissues, including use in vascular and ophthalmic procedures. It is not indicated to be used in the central circulatory system and central nervous system.

ACTIONS

Nylon Monofilament Suture causes a minimal acute inflammatory reaction in tissue, which is followed by a gradual encapsulation of the suture by fibrous connective tissue. While polyamide is not absorbed, progressive hydrolysis of the polyamide *in vivo* may result in gradual loss of tensile strength over time.

CONTRAINDICATION

Due to the gradual loss of tensile strength which may occur over prolonged periods *in vivo*, the nylon sutures should not be used where permanent retention of tensile strength is required.

This device is contraindicated to persons that could have any allergic reaction to the medical device material.

WARNINGS

Users must be acquainted with surgical procedures and practices involving non absorbable sutures before using the Nylon Monofilament Suture to close incisions for the risk of wound dehiscence may vary according to the site of application and the suture material used. As which any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculi formation. The suture should only be used by or under surgeon's orientation according to applicable surgical techniques and it should only be manipulated by qualified people. Do not use the suture for other purposes, just the one indicated by the manufacturer. When handling this suture be sure of doing it aseptically. Acceptable surgical practices must be followed for management of infected or contaminated wounds. Do not re-sterilize. Discard the open packages and unused sutures. Discard left over. Do not use if packing is violated or damp. Do not use the product beyond the validity period. Disposal of used sutures: Risk of biological contamination - adopt internal disposal procedures, taking into account the risk of biological contamination. Disposal of used and non-used sutures: Perforating and cutting materials - adopt internal disposal procedures for perforating and cutting materials.

PRECAUTIONS

In handling this or any other surgical suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Broken needles may result in extended or additional surgeries or residual foreign bodies. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in recipients for potentially contaminated material.

USE

Use the necessary amount of suture in agreement with the extension of the surgical wound and the technique adopted by the surgeon. The size and diameter of suture and needle used for each type of surgical intervention are defined as per medical criteria.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include, gradual loss of tensile strength over time wound dehiscence, infection, acute inflammatory tissue reaction, transitory local irritation at the wound site and calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs.

The surgeon must inform the patient to look for medical assistance in case of any suspect of adverse reactions symptoms.

STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30°C (59°F to 86°F), protected from direct sunlight, heat sources, dust and humidity.

HOW SUPPLIED

The Nylon Monofilament Suture is supplied dyed in blue, black or undyed, in sterile envelope, in U.S.P sizes 12-0 through 10 (metric sizes 0.01 through 12). Suture is supplied in several lengths, without needle or with 1 or 2 non traumatic needles, according to each type of surgical intervention in boxes with 12, 24 or 36 units.

EXPIRATION

Five years from manufacturing date. See box and envelope.

PACKAGES WITH 12 (TWELVE), 24 (TWENTY FOUR) OR

36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.

SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES

DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN - DESTROY AFTER USING

STERILIZATION METHOD: GAMMA RADIATION OR ETHYLENE OXIDE

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CNPJ 02.370.649/0001-20 - Fábrica: Tel.: 55 19 3646-7300

Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

REGISTRO ANVISA N.º:

NYLON MONOFILAMENT NO NEEDLE: 10372360014

NYLON MONOFILAMENT WITH NEEDLE: 10372360009