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Responsable Técnico: MARINA GIL CORSI - CRF/SP 94.142

REGISTRO POLIGLACTINA CON AGUJA NA ANVISA Nº: 10372360024

INSTRUCTIONS FOR USE OF POLIGLACTIN SUTURE (PGL)

SYNTHETIC ABSORBABLE SURGICAL SUTURE

DESCRIPTION

Polyglactin Suture (PGL) is a synthetic, absorbable, sterile, multifilament surgical suture, composed of a copolymer made from 90% Glicolide and 10% L-Lactide, and is coated with a mixture composed of equal parts of copolymer of Polyglactin 370 (30% Glicolide and 70% L-Lactide) with Calcium Stearate. The copolymers used in this product have been found to be nonantigenic, nonpyrogenic and elicit only a mild tissue reaction during absorption. Polyglactin Suture is available dyed violet to enhance its visibility on the tissue or on its natural color (beige). The pigment used in the violet suture is D&C #2 approved by the FDA. This product meets specifications set forth for Absorbable Surgical Sutures by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia (EP) and for the Ethylene Oxide NBR 13904 - Filaments for Surgical Sutures. Polyglactin Suture (PGL) is sterilized by Ethylene Oxide. Polyglactin Suture (PGL) meet all requirements established by the norms for absorbable surgical sutures, except the following diameter:

MAXIMUM SUTURE OVERSIZE IN DIAMETER (mm) U.S.P	
U.S.P SUTURE SIZE DESIGNATION	MAXIMUM OVERSIZE "mm"
6-0	0,008
5-0	0,014
4-0	0,017
3-0	0,018
2-0	0,004
0 and 1	0,020

INDICATION

Polyglactin Suture (PGL) is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not indicated for use in cardiovascular or neurological tissues.

ACTIONS

Polyglactin Suture (PGL) elicits a minimal acute inflammatory reaction in tissue followed by a gradual encapsulation of the suture by the fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of Suture PGL occurs by means of hydrolysis, where the polymer degrades to glycolic and lactic acids which is subsequently absorbed and metabolized in the body. Absorption begins as loss of tensile strength followed by a loss of mass. Absorption studies *in vitro* indicate that after 14 days the tensile strength is approximately 60% of its initial tensile strength. The suture absorption is completed between 56-70 days after the implant.

CONTRAINDICATION

Polyglactin Suture (PGL) being absorbable should not be used where extended approximation of tissue under stress is required.

WARNINGS

Users must be acquainted with surgical procedures and practices involving absorbable sutures before using Polyglactin Suture (PGL) to close incisions for the risk of wound dehiscence may vary according to the site of application and the suture material used. The use of this suture may be unsuitable for elder, undemourished, weak or other patients whose conditions may slow down the wound from skinning over. As this is an absorbable suture material the use of supplemental non absorbable sutures should be considered by the surgeon in the closure of wounds which may undergo expansion, stretching or distension, or which may require additional support. Polyglactin Suture as an absorbable suture it may act transiently as a foreign body. As with any foreign body prolonged contact of any suture with saline solution, such as those found at the biliary and urinary canals, may result in calculi formation. The suture should only be used for or under surgeon's orientation according to applicable surgical techniques and it should only be manipulated by qualified people. Do not use the suture for other purposes, just the one indicated by the manufacturer. When handling this suture be sure of doing it aseptically. Acceptable surgical practices must be followed to treat any infected or contaminated wounds. Do not resterilize. Discard the open packages and unused sutures. Discard left over. Do not use if packing is violated or damp. Do not use the product beyond the validity period. Disposal of used sutures: Risk of biological contamination – adopt internal disposal procedures, taking into account the risk of biological contamination. Disposal of used and non-used sutures: Perforating and cutting materials – adopt internal disposal procedures for perforating and cutting materials.

PRECAUTIONS

Under certain circumstances, especially in orthopedic procedures, the immobilization of joints by external support may be employed, if the surgeon finds it necessary. It is necessary to reflect on the use of absorbable sutures in tissues with low supply of blood, for it may occur extrusion of the suture and the slowing down of the absorption. In the handling this or any other suture material care should be taken to avoid damage from handling. Avoid the crimping and the crushing of the needles by the inadequate handling of tweezers and needle holders. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens. Dispose used needles in recipients for potentially contaminated material.

USE

Use the necessary amount of suture according to the wound extension and the technique adopted by the surgeon. The size and diameter of suture and needle used for each type of surgical intervention are defined as per medical criteria.

ADVERSE REACTIONS

Adverse reactions connected with the use of this absorbable synthetic suture include: wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, tissue minimal acute inflammatory reaction, localized irritation when skin suture are left in place for greater than 7 days, suture extrusion and slow absorption in tissues with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation the wound site.

STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30° C, protected from direct sunlight, heat sources, humidity and dust.

HOW SUPPLIED

Polyglactin Suture (PGL) is supplied in violet and in its natural color (beige); in U.S.P. sizes 10-0 through 2 (metric sizes 0,2 through 5) Suture is supplied in a variety of lengths, attached to many types of non traumatic needles, according to each type of surgical intervention; in boxes with 12, 24 or 36 units of sterile pouches.

EXPIRATION

Five years from manufacturing date. See box and envelope.

PACKAGES WITH 12 (TWELVE), 24 (TWENTY FOUR) OR

36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.

SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES

DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN - DESTROY AFTER USING

STERILIZATION METHOD : ETHYLENE OXIDE

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CNPJ 02.370.649/0001-20 - Fábrica: Tel.: 55 19 3646-7300

Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

ANVISA REGISTRATION POLYGLACTIN SUTURE WITH NEEDLE: Nº 10372360024