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INSTRUCTIONS FOR USE OF POLYPROPYLENE SUTURE

NON ABSORBABLE SURGICAL SUTURE

DESCRIPTION

Polypropylene Suture is a synthetic, non absorbable, sterile, monofilament suture, composed of isostatic polypropylene polymer of high molecular weight. The suture is presented colored blue to enhance its visibility on the tissue. The pigments used are D&C blue nbr 6, solvent blue nbr 104 or CU Phtalocyanine blue, approved by the FDA. The Polypropylene Suture meets such specifications as set forth for Non Absorbable Surgical Sutures by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia (EP) and for the ABNT norm NBR 13904 - Filaments for Surgical Sutures. The Polypropylene Suture is sterilized by Ethylene Oxide. Polypropylene Suture meet all requirements established by the norms for non absorbible surgical sutures, except the following diameter:

MÁXIMUM SUTURE OVERSIZE IN DIAMETER (mm)	
U.S.P SUTURE SIZE DESIGNATION	MAXIMUM OVERSIZE "mm"
7-0	0.007

INDICATION

Polypropylene Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

ACTIONS

Polypropylene Suture elicits a minimal acute inflammatory reaction in tissue which is followed by gradual encapsulation of the suture by fibrous connective tissue. Polypropylene Suture is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. As a monofilament, Polypropylene Suture resists involvement in infection and has been successfully employed in contaminated and infected wounds to eliminate or minimize later sinus formation and suture extrusion. The lack of adherence to tissues has facilitated the use of Polypropylene Suture as pull-out suture.

CONTRAINDICATION

This device is contraindicated to persons that could have any allergic reaction to the medical device material.

WARNINGS

Users should be familiar with surgical procedures and techniques involving non absorbable sutures before employing Polypropylene Suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. As which any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculi formation. of infected or contaminated wounds. Acceptable surgical practice should be followed for the management under surgeon's orientation according to applicable surgical techniques and it should only be manipulated by qualified people. Do not use the suture for other purposes, just the one indicated by the manufacturer. When handling this suture be sure of doing it aseptically. Acceptable surgical practice should be followed for the management of infected or contaminated wounds. Do not re-sterilize. Discard the open and unused packages. Discard left over. Do not use if packing is violated or damp. Do not use the product beyond the validity period. Disposal of used sutures: Risk of biological contamination - adopt internal disposal procedures, taking into account the risk of biological contamination. Disposal of used and non-used sutures: Perforating and cutting materials - adopt internal disposal procedures for perforating and cutting materials.

PRECAUTIONS

In handling this or any other surgical suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Braking of needles may cause surgeries to take longer and additional surgeries and foreign matters left inside the body. Inadvertent piercing of users with contaminated needles may cause transmission of blood borne pathogens. Users must be careful when handling surgical needles to avoid undesirable perforations. Dispose used needles in recipients for potentially contaminated material.

USE

Use the necessary amount of suture in agreement with the extension of the surgical wound and the technique adopted by the surgeon. The size and diameter of suture and the needle used for each type of surgical intervention are defined as per medical criteria.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, infection, minimal acute inflammatory tissue reaction, transitory local irritation at the wound site, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs. The surgeon must inform the patient to look for medical assistance in case of any suspect of adverse reactions symptoms.

STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30°C (59°F to 86°F), protected from direct sunlight, heat sources, humidity and dust.

HOW SUPPLIED

Polypropylene Suture is supplied in blue color, in sterile envelopes, in U.S.P sizes 10-0 through 6 (metric sizes 0.2 through 8). Suture is supplied in several lengths, attached with 1 or 2 non traumatic needles, attached or not to TFE pledgets, according to each type of surgical intervention; in boxes with 12, 24 or 36 units.

EXPIRATION

Five years from manufacturing date. See box and envelope.

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PACKAGES WITH 12 (TWELVE), 24 (TWENTY FOUR) OR
36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.
SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES
DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN - DESTROY AFTER USING
STERILIZATION METHOD: ETHYLENE OXIDE

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CNPJ 02.370.649/0001-20 - Fábrica: Tel.: 55 19 3646-7300

Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

REGISTRO POLIPROPILENO CON AGUJA ANVISA Nº: 10372360006

REGISTER POLYPROPYLENE WITH NEEDLE ANVISA Nº: 10372360006