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INSTRUCTIONS FOR USE OF SILK SUTURE

NON ABSORBABLE SURGICAL SUTURE

DESCRIPTION

Silk Suture is a non absorbable, sterile, multifilament, braided or twisted suture, made of an organic protein called fibroin. This protein is derived from the domesticated species Bombyx Mori (B. Mori) of the Bombycidae family. The suture is dyed black or blue to enhance its visibility in the tissue. The coloring pigments are approved by the FDA being Logwood Black or Sulphur Black 1, for black silk and the Methylene Blue for blue silk. The home and foreign suppliers may use other pigments since they are biologically compatible, approved by the FDA and not having toxicity. This suture is coated with silicone. This suture meets the specifications set forth for Non Absorbable Surgical Sutures by the Brazilian Pharmacopoeia, by the United States Pharmacopoeia (USP), by the European Pharmacopoeia (EP) and for the ABNT norm NBR 13904 - Filaments for Surgical Sutures. It is sterilized by Gamma Radiation (Cobalt 60) or Ethylene Oxide.

INDICATION

Silk Suture is indicated for use in general soft tissue approximation and/or ligation, including use in vascular, ophthalmic and neurological procedures. It is not indicated to be used in the central circulatory system and central nervous system. The Silk Suture is not indicated to be used in the neurological areas: brain, meninges and spinal cord. The neurological indication of Silk Suture is other than Brain, meninges and spinal cord, like the nervous ramifications in the members.

ACTIONS

Silk Suture elicits an acute inflammatory reaction in tissue which is followed by gradual encapsulation of the suture by the fibrous connective tissue. While the Silk Suture is not absorbed, progressive degradation of the proteinaceous silk fiber *in vivo* may result in gradual loss of all the suture tensile strength over time.

CONTRA INDICATION

Silk Suture is contraindicated in patients with known sensitivities or allergies to silk. Due to gradual loss of tensile strength which may occur over prolonged periods *in vivo* Silk Suture should not be used where permanent retention of tensile strength is required.

This device is contraindicated to persons that could have any allergic reaction to the medical device material.

WARNINGS

Users should be familiar with surgical procedures and techniques involving non absorbable sutures before employing Silk Suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculi formation. The suture should be used only by or under surgeon's orientation according to applicable surgical techniques and it should only be manipulated by qualified people. Acceptable surgical practice should be followed. Do not use the suture for other purposes, just the one indicated by the manufacturer. When handling this suture be sure of doing it aseptically. Acceptable surgical practice must be followed for the management of infected or contaminated wounds. Do not re-sterilize. Discard open and unused packages. Discard left over. Do not use if packing is violated or damp. Do not use the product beyond the validity period. Disposal of used sutures: Risk of biological contamination - adopt internal disposal procedures, taking into account the risk of biological contamination. Disposal of used and non-used sutures: Perforating and cutting materials - adopt internal disposal procedures for perforating and cutting materials.

PRECAUTIONS

In handling this or any other surgical suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Broken needles may cause longer and additional surgeries and foreign matters left inside the body. Inadvertent piercing of users with contaminated needles may cause transmission of blood borne pathogens. Dispose used needles in recipients for potentially contaminated material.

USE

Use the necessary amount of suture in agreement with the extension of the surgical wound and the technique adopted by the surgeon. The size and diameter of suture and the needle used for each type of surgical intervention are defined as per medical criteria.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, gradual loss of tensile strength over time, allergic response in patients sensitive to silk, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, transitory local irritation at the wound site.

The surgeon must inform the patient to look for medical assistance in case of any suspect of adverse reactions symptoms.

STORAGE AND TRANSPORTATION

Transport and store the product in its original packing, at temperatures from 15°C to 30°C (59°F to 86°F), protected from direct sunlight, heat sources, humidity and dust.

HOW SUPPLIED

Silk Suture is supplied in blue or black colors, in sterile envelopes, in U.S.P sizes 10-0 through 6 (metric sizes 0.2 through 8). Suture is supplied in several lengths, without or with 1 or 2 non traumatic needles, according to each type of surgical intervention; in boxes with 12, 24 or 36 units.

EXPIRATION

Five years from manufacturing date. See box and envelope.

PACKAGES WITH 12 (TWELVE), 24 (TWENTY FOUR) OR

36 (THIRTY SIX) UNITS OF STERILE ENVELOPES.

SINGLE USE PRODUCT - DISCARD ANY OPEN AND NOT USED TAPES

DO NOT RE-STERILIZE - STERILE - REUSE IS FORBIDDEN - DESTROY AFTER USING

STERILIZATION METHOD: GAMMA RADIATION OR ETHYLENE OXIDE

BRASUTURE IND. COM. IMP. EXP. LTDA.

TRADE MARK / MADE IN BRAZIL

Rua Vereador José Vasconcellos dos Reis, n.º 642 - Distrito Industrial
CEP 13790-000 - São Sebastião da Gramma - SP - Brasil
CNPJ 02.370.649/0001-20 - Fábrica: Tel.: 55 19 3646-7300
Technical Responsible: MARINA GIL CORSI - CRF/SP 94.142

REGISTRO ANVISA N.º:

SILK SUTURE NO NEEDLE: 10372360016

SILK SUTURE WITH NEEDLE: 10372360011